

Automated Haematology Analyser  
**XN-L Series**

Leading Haematology  
for Better Patient Care



# Upgrade your 3PD to 6PD with Sysmex

## Why 5PD?

### 3PD Analyser

Neutrophils

Mixed

Lymphocytes

### 5PD Analyser

Neutrophils

Eosinophils

Basophils


Monocytes

Lymphocytes

Five part differential analysers count monocytes, eosinophils and basophils separately rather than as a mixed population in three part differential analysers. These individual WBC parameters reveal more valuable information to clinicians to support diagnostic and treatment decisions.


## Why Sysmex 6PD?

### 3PD Analyser



SFL

### Sysmex 6PD Analyser



SSC

added-value  
IG<sup>1</sup>

Examples of WBC Flags:  
**WBC Abn Scattergram**  
**Neutropenia**  
**Neutrophilia**  
**Lymphopenia**  
**Lymphocytosis**  
**Monocytosis**  
**Eosinophilia**  
**Leukocytopenia**  
**Leukocytosis**  
**Blasts/Abn Lympho?**  
**Left Shift?**  
**Atypical Lympho?**  
**NRBC?**

## Benefits of Sysmex 6PD Analyser

### 1) Assessment is not solely based on cell size

Unlike three part differential and some five part differential analysers, fluorescent flow cytometry measures not only the cell size, but also the intracellular information and nucleic acid content. This produces a highly accurate differential count in EDTA blood samples even as the cell size changes during normal storage.

### 2) Identification of immature cells

Identification of immature cells is possible with the XN-L as immature cells have a higher nucleic acid content. This has made the generation of six part differential, the IG count possible. The precision of an automated IG count increases laboratory efficiency by reducing manual counts.

### 3) Superior flagging system for abnormal cells

As compared to a flagging system based entirely on cell size on a three part differential analyser, the XN-L provides more detailed and specific flagging for abnormal cells.

# Introducing XN-L Series

XN-L Series is the latest compact fully-automated 6-part differential haematology analyser from Sysmex.

It is designed to meet today's laboratory needs by providing enhanced clinical values that only high-end models were previously able to provide; delivering improved operational efficiency in the laboratories.

XN-L Series is available in four models with different aspiration modes.



## XN-550

Sampler that allows both open and closed tube analysis



## XN-450

Closed tube analysis



## XN-330 / XN-350

Open tube analysis

# A completely new compact 6-part differential haematology analyser that caters for what matters most in your lab.

## 1 Ease of use & peace of mind

Multi-coloured touch screen display embedded with intuitive graphical user interface to support your routine tasks.




A reagent management system that calculates remaining tests and expiry dates to allow more efficient tracking of reagent usage.



## 2 Space saving

Compact design with an in-built IPU that fits easily to any laboratory bench or table with a smaller footprint.



## 3 Flexibility

Only **25µL** of aspiration volume

Requires only 25 µl of aspiration volume for whole blood mode and low WBC mode.

### Auto-diluent dispensing function

Auto-diluent dispensing function is available for dilution of samples

## 4 Improved workflow and TAT

Up to **70** samples/hour\*

Better TAT with one of the highest throughput compact 6PD analysers in the market.

Sampler available for a truly walk-away system\*\*




Left drawer



Right drawer

Sampler comes with two drawers to allow continuous loading of samples

Integrated auto Repeat, Rerun and Reflex\*\* of sample




\*Only available on XN-350 / XN-450 / XN-550

\*\*Only available on XN-550

# More than just CBC + DIFF

## Core Technology of XN-L Series

XN-L Series utilises the same core technology as high-end haematology analysers - the proven and reliable laser flow cytometry - to count and analyse cells. After the cells are irradiated by the laser beam, the forward scattered light (FSC), side scattered light (SSC) and side fluorescent light (SFL) of the cells are analysed. The three signals are used to differentiate and count cells with the help of unique digital technology and algorithms.



## Standard application

added-value  
**IG<sup>1,2</sup>**

**WDF Scattergram**

**Added value:**  
A true 6-part differential count (IG):

- Immature granulocytes (IG) with every differential analysis aids in the early prediction of infection and inflammation.
- **No additional reagents required.**

## Optional applications\*


added-value  
**Low  
WBC**

### Added value:

Analysis of low WBC samples in L-WBC mode:

- Samples with low white blood cell count can be measured accurately by counting twice as many cells; providing reliable results that aid in chemotherapy monitoring.


added-value  
**Body  
Fluid**<sup>3,4</sup>



### Added value:

Fully-automated body fluid analysis in BF mode:


- 2-part differential body fluid analysis includes MN (mononuclear) and PMN (polymorphonuclear) cell population to aid in the distinction between viral and bacterial infection.



- **No additional reagents required.**
- No special sample preparation required.

added-value  
**RET<sup>5</sup>**

### RET Scattergram




### Added value:

Indices of erythropoiesis (RET, RET-He, IRF):

- RET-He (reticulocytes haemoglobin) and IRF (immature reticulocytes fraction) aids in monitoring of RBC production.
- RET-He (reticulocytes haemoglobin) aids in differentiation between functional and classical iron deficiency and monitoring of EPO and/or IV iron therapy.

added-value  
**IPF<sup>6,7,8</sup>**

### PLT-O Scattergram



### Added value:

Indices of thrombopoiesis (IPF):


- IPF (Immature platelet fraction) aids in differential diagnosis of thrombocytopenic disorders and is an early predictor of platelet recovery.

\*Only available on XN-350 / XN-450 / XN-550

# What's more

## Extended Connectivity

- XN-L supports extended connectivity to meet a laboratory's various connectivity requirements.



## Standardisation across Sysmex platforms

- Compatibility with the XN-Series through common reagents, controls and technology allows for consistency across laboratory operations.

### 1 Common reagents and QC materials



Reagents and QC materials utilised with the XN-Series can be used on the XN-L Series.

### 2 Uniform user interface



The XN-L analysers use the same user interface as the XN-Series.

### 3 Identical core technology and clinical parameters



Enhanced clinical parameters makes the XN-L Series the perfect backup for 24/7 laboratories.

# Specifications



XN-330



XN-350



XN-450



XN-550

|  |  |   |   |   |
|--|--|---|---|---|
| <b>Principles &amp; Technologies</b>       | <p><b>Fluorescent Flow Cytometry method:</b><br/>WBC, DIFF<br/>RET, IPF, 2-part differential for body fluid analysis (except XN-330)</p> <p><b>Hydrodynamic focusing DC detection method:</b><br/>PLT-I (Impedance), RBC, HCT</p> <p><b>Cyanide-free SLS-haemoglobin method:</b><br/>HGB</p> |   |   |   |
| <b>Parameters</b>                          | <p><b>28 Parameters</b><br/><u>Whole blood / Pre-dilution mode:</u><br/>WBC, RBC, HGB, HCT, MCV, MCH, MCHC, PLT, RDW-SD, RDW-CV, PDW, MPV, P-LCR, PCT, NEUT#, LYMPH#, MONO#, EO#, BASO#, NEUT%, LYMPH%, MONO%, EO%, BASO%, IG#, IG%, MicroR, MacroR</p>                                      |   |   |   |
| <b>Optional Parameters (Except XN-330)</b> | <p><b>14 RET &amp; IPF Parameters</b><br/>RET#, RET%, IRF, LFR, MFR, HFR, RET-He, PLT-O, IPF#, IPF, RBC-He, Delta-He, HYPO-He, HYPER-He</p> <p><b>7 Body Fluids Parameters</b><br/><u>Body Fluid mode</u><br/>WBC-BF, RBC-BF, MN#, PMN#, MN%, PMN%, TC-BF#</p>                               |   |   |   |
| <b>Throughput</b>                          | <p>CBC :<br/>up to 60 samples/hour</p> <p>CBC + DIFF:<br/>up to 60 samples/hour</p>  | <p>CBC: up to 70 samples/hour<br/>CBC + DIFF: up to 70 samples/hour<br/>CBC + DIFF + RET: up to 35 samples/hour</p> <p>Body fluid mode: up to 30 samples/hour</p> |   |   |
| <b>Sample Aspiration Volume</b>            | <p>Whole blood mode: 25 µl<br/>Pre-dilution mode: 70 µl</p>  | <p>Whole blood / Low WBC mode: 25 µl<br/>Pre-dilution mode: 70 µl<br/>Body fluid mode: 70 µl</p>  |   |   |
| <b>Data Storage</b>                        | <p>Results: 100,000 samples<br/>Patient information: 10,000 records<br/>QC files: 99 files / analyser<br/>QC plots : 300 plots / file<br/>Reagent replacement history: 5,000 records<br/>Maintenance history: 5,000 records</p>  |   |   |   |
| <b>Dimensions and Weight</b>               | <p>Width : 450 mm<br/>Depth : 460 mm<br/>Height : 510 mm<br/>Approx. 35 kg</p>   | <p>Width: 450 mm<br/>Depth: 460 mm<br/>Height: 510 mm<br/>Approx. 35 kg</p>   | <p>Width: 450 mm<br/>Depth: 460 mm<br/>Height: 440 mm<br/>Approx. 35 kg</p> | <p>Width: 450 mm<br/>Depth: 660 mm<br/>Height: 450 mm<br/>Approx. 53 kg</p> <p><b>Monitor</b><br/>Width: 267 mm<br/>Depth: 205 mm<br/>Height: 240 mm<br/>Approx. 3 kg</p> |

## References :-

1. Sang Ook Ha, Sang Hyuk Park, So Hee Park, Jae Seok Park, Jin Won Huh, Chae-Man Lim, Younsuck Koh, Sang-Bum Hong & Seongsoo Jang (2015) Fraction of immature granulocytes reflects severity but not mortality in sepsis, *Scand J Clin Lab Invest*, 75:1, 36-43, DOI: 10.3109/00365513.2014.965736
2. Nierhaus A, Klatte S, Linssen J, Eismann NM, Wichmann D, Hedke J, et al. Revisiting the white blood cell count: immature granulocytes count as a diagnostic marker to discriminate between SIRS and sepsis - a prospective, observational study. *BMC Immunol* 2013;14(1):8.
3. Fleming C, Russcher H, Lindemans J, De Jonge R. Clinical relevance and contemporary methods for counting blood cells in body fluids suspected of inflammatory disease. Vol. 53, *Clinical Chemistry and Laboratory Medicine*. 2015. p. 1689-706.
4. Seghezzi M, Buoro S, Manenti B, Mecca T, Ferrari R, Zappalà G, et al. Optimization of Cellular analysis of Synovial Fluids by optical microscopy and automated count using the Sysmex XN Body Fluid Mode. *Clin Chim Acta*. 2016;462:41-8.
5. Weimann A, Cremer M, Hernáiz-Driever P, Zimmermann M. Delta-He, Ret-He and a new diagnostic plot for differential diagnosis and therapy monitoring of patients suffering from various disease-specific types of anemia. *Clin Lab*. 2016;62(4):667-77.
6. Dadu T, Sehgal K, Joshi M, Khodaiji S. Evaluation of the immature platelet fraction as an indicator of platelet recovery in dengue patients. *Int J Lab Hematol*. 2014;36(5):499-504.
7. Sakuragi M, Hayashi S, Maruyama M, Kabutomori O, Kiyokawa T, Nagamine K, et al. Clinical significance of IPF% or RP% measurement in distinguishing primary immune thrombocytopenia from aplastic thrombocytopenic disorders. *Int J Hematol*. 2015;101(4):369-75.
8. van der Linden N, Klinkenberg LJJ, Meex SJR, Beckers EAM, de Wit NCJ, Prinzen L. Immature platelet fraction measured on the Sysmex XN hemocytometer predicts thrombopoietic recovery after autologous stem cell transplantation. *Eur J Haematol*. 2014;93(2):150-6.



---

**Sysmex Asia Pacific Pte Ltd**  
Tel +65 6221-3629 Fax +65 6221-3687  
[www.sysmex-ap.com](http://www.sysmex-ap.com)

**Sysmex India Pvt. Ltd**  
Tel +91 (22) 6112-6666  
[www.sysmex.co.in](http://www.sysmex.co.in)

**PT Sysmex Indonesia**  
Tel +62 (21) 3002-6688 Fax +62 (21) 3002-6699  
[www.sysmex.co.id](http://www.sysmex.co.id)

**Sysmex (Malaysia) Sdn Bhd**  
Tel +60 (3) 5637-1788 Fax +60 (3) 5637-1688  
[www.sysmex.com.my](http://www.sysmex.com.my)

**Sysmex New Zealand Ltd**  
Tel +64 (9) 630-3554  
[www.sysmex.co.nz](http://www.sysmex.co.nz)

**Sysmex (Thailand) Co., Ltd**  
Tel +66 (2) 032-2536 Fax +66 (2) 116-5396  
[www.sysmex.co.th](http://www.sysmex.co.th)

**Sysmex Philippines Inc.**  
Tel +63 (2) 621-2460 Fax +63 (2) 621-2432  
[www.sysmex.com.ph](http://www.sysmex.com.ph)

**Sysmex Vietnam Co., Ltd**  
Tel +84 (8) 3997-9400 Fax +84 (8) 3997-9405  
[www.sysmex.com.vn](http://www.sysmex.com.vn)

**Sysmex Vietnam Co., Ltd (Hanoi Branch)**  
Tel +84 (4) 3776-7020 Fax +84 (4) 3776-7022  
[www.sysmex.com.vn](http://www.sysmex.com.vn)

**Sysmex Australia Pty Ltd**  
Tel +61 (2) 9016-3040  
[www.sysmex.com.au](http://www.sysmex.com.au)